

兔抗 CORO1C 多克隆抗体

中文名称：兔抗 CORO1C 多克隆抗体

英文名称： Anti-CORO1C rabbit polyclonal antibody

别名： HCRNN4

相关类别： 一抗

储存： 冷冻（-20℃）

抗原： CORO1C

宿主： Rabbit

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene encodes a member of the WD repeat protein family. WD repeats are minimally conserved regions of approximately 40 amino acids typically bracketed by gly-his and trp-asp (GH-WD), which may facilitate formation of heterotrimeric or multiprotein complexes. Members of this family are involved in a variety of cellular processes, including cell cycle progression, signal transduction, apoptosis, and gene regulation. Three transcript variants encoding two different isoforms have been found for this gene.
Applications:	ELISA, WB, IHC

Name of antibody:	CORO1C
Immunogen:	Fusion protein of human CORO1C
Full name:	coronin, actin binding protein, 1C
Synonyms :	HCRNN4
SwissProt:	Q9ULV4
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human liver cancer and human colon cancer
IHC Recommend dilution:	100-300
WB Predicted band size:	53 kDa
WB Positive control:	Lovo cells
WB Recommended dilution:	500-2000



