

## 兔抗 CORO6 多克隆抗体

中文名称：兔抗 CORO6 多克隆抗体

英文名称：Anti-CORO6 rabbit polyclonal antibody

别名：coronin 6

相关类别：一抗

储存：冷冻（-20℃）

抗原：CORO6

宿主：Rabbit

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	Coronin-6 also known as coronin-like protein E (Clipin-E) is a protein that in humans is encoded by the CORO6 gene. Coronin-6 belongs to coronin family which is an actin binding protein. Human CORO6 gene is located on chromosome 17 on the cytogenetic band 17 p11.2. Gene CORO6 is well conserved across domain of eukaryotic organisms from animal to fungi.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	CORO6
<b>Immunogen:</b>	Full length fusion protein
<b>Full name:</b>	coronin 6

<b>SwissProt:</b>	Q6QEF8
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human thyroid cancer
<b>IHC Recommend dilution:</b>	25-100

