

CNOT2 抗原（重组蛋白）

中文名称： CNOT2 抗原（重组蛋白）

英文名称： CNOT2 Antigen (Recombinant Protein)

别名： CCR4-NOT transcription complex subunit 2; NOT2; CDC36; NOT2H; HSPC131

储存： 冷冻（-20℃）

相关类别： 抗原

概述

Fusion protein corresponding to a region derived from 341-540 amino acids of human CNOT2

技术规格

Full name:	CCR4-NOT transcription complex subunit 2
Synonyms:	NOT2; CDC36; NOT2H; HSPC131
Swissprot:	Q9NZN8
Gene Accession:	BC002597
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene encodes a subunit of the multi-component CCR4-NOT complex. The CCR4-NOT complex regulates mRNA synthesis and degradation and is also thought to be involved in mRNA splicing, transport and localization. The encoded protein interacts with histone deacetylases and functions as a repressor of polymerase II transcription. Alternatively spliced transcript variants have been observed for this gene.

mlbio 梅联生物
Good elisakit producers