

兔抗 CORO7 多克隆抗体

中文名称：兔抗 CORO7 多克隆抗体

英文名称： Anti-CORO7 rabbit polyclonal antibody

别名： CRN7; POD1; 0610011B16Rik

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： CORO7

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	Coronins, such as CORO7, constitute an evolutionarily conserved family of WD-repeat actin-binding proteins. CORO7 plays a role in Golgi complex morphology and function.
Applications:	ELISA, IHC
Name of antibody:	CORO7
Immunogen:	Fusion protein corresponding to a region derived from 726-925 amino acids of human CORO7
Full name:	coronin 7
Synonyms:	CRN7; POD1; 0610011B16Rik
SwissProt:	P57737

ELISA Recommended dilution:	5000-10000
IHC Positive control:	Human liver cancer;Human esophagus cancer
IHC Recommended dilution:	25-100



