

CNPY2 抗原（重组蛋白）

中文名称：CNPY2 抗原（重组蛋白）

英文名称：CNPY2 Antigen (Recombinant Protein)

别名：MSAP; TMEM4; ZSIG9; HP10390

储存：冷冻（-20℃）

相关类别：抗原

概述

Full length fusion protein

技术规格

Full name:	canopy FGF signaling regulator 2
Synonyms:	MSAP; TMEM4; ZSIG9; HP10390
Swissprot:	Q9Y2B0
Gene Accession:	BC065015
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	Positive regulator of neurite outgrowth by stabilizing myosin regulatory light chain (MRLC). It prevents MIR-mediated MRLC ubiquitination and its subsequent proteasomal degradation.