

兔抗 ASF1A 多克隆抗体

中文名称: 兔抗 ASF1A 多克隆抗体

英文名称: Anti-ASF1A rabbit polyclonal antibody

别名: CIA; CGI-98; HSPC146

储存: 冷冻 (-20℃)

抗原: ASF1A

宿主: Rabbit

反应种属: Human, Mouse

相关类别: 一抗

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	This gene encodes a member of the H3/H4 family of histone chaperone proteins and is similar to the anti-silencing function-1 gene in yeast. The protein is a key component of a histone donor complex that functions in nucleosome assembly. It interacts with histones H3 and H4, and functions together with a chromatin assembly factor during DNA replication and repair.
Applications:	ELISA, WB
Name of antibody:	ASF1A
Immunogen:	Synthetic peptide of human ASF1A
Full name:	anti-silencing function 1A histone chaperone

Synonyms:	CIA; CGI-98; HSPC146
SwissProt:	Q9Y294
ELISA Recommended dilution:	2000-5000
WB Predicted band size:	23 kDa
WB Positive control:	HEPG2 and PC-3 cell lysates
WB Recommended dilution:	200-1000

