

## 兔抗 ASH2L 多克隆抗体

中文名称：兔抗 ASH2L 多克隆抗体

英文名称：Anti-ASH2L rabbit polyclonal antibody

别名：ASH2; Bre2; ASH2L1; ASH2L2

相关类别：一抗

储存：冷冻（-20℃）

抗原：ASH2L

宿主：Rabbit

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	Component of the Set1/Ash2 histone methyltransferase (HMT) complex, a complex that specifically methylates 'Lys-4' of histone H3, but not if the neighboring 'Lys-9' residue is already methylated. As part of the MLL1/MLL complex it is involved in methylation and dimethylation at 'Lys-4' of histone H3. May function as a transcriptional regulator. May play a role in hematopoiesis.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	ASH2L
<b>Immunogen:</b>	Fusion protein of human ASH2L

<b>Full name:</b>	ASH2 like, histone lysine methyltransferase complex subunit
<b>Synonyms:</b>	ASH2; Bre2; ASH2L1; ASH2L2
<b>SwissProt:</b>	Q9UBL3
<b>ELISA Recommended dilution:</b>	1000-5000
<b>IHC positive control:</b>	Human colon cancer and Human thyroid cancer
<b>IHC Recommend dilution:</b>	50-200

