

COMMD6 抗原（重组蛋白）

中文名称：COMMD6 抗原（重组蛋白）

英文名称：COMMD6 Antigen (Recombinant Protein)

别名：Acrg

储存：冷冻（-20℃）

相关类别：抗原

概述

Full length fusion protein

技术规格

Full name:	COMM domain containing 6
Synonyms:	Acrg
Swissprot:	Q7Z4G1
Gene Accession:	BC117391
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	COMMD6 belongs to a family of NF-kappa-B (see RELA; MIM 164014)-inhibiting proteins characterized by the presence of a COMM domain (see COMMD1; MIM 607238) (de Bie et al., 2006 [PubMed 16573520])