

兔抗 ASMTL 多克隆抗体

中文名称：兔抗 ASMTL 多克隆抗体

英文名称：Anti-ASMTL rabbit polyclonal antibody

别名：acetylserotonin O-methyltransferase like; ASMTL; ASMTLX; ASMTLY

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：ASMTL

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	The protein encoded by this gene has an N-terminus that is similar to the multicopy associated filamentation (maf) protein of <i>Bacillus subtilis</i> and to orfE of <i>Escherichia coli</i> , while the C-terminus is similar to N-acetylserotonin O-methyltransferase. This gene is located in the pseudoautosomal region 1 (PAR1) of X and Y chromosomes. Three transcript variants encoding different isoforms have been found for this gene.
Applications:	ELISA, IHC
Name of antibody:	ASMTL
Immunogen:	Fusion protein of human ASMTL

Full name:	acetylserotonin O-methyltransferase like
Synonyms:	ASTML; ASMTLX; ASMTLY
SwissProt:	O95671
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human tonsil and Human thyroid cancer
IHC Recommend dilution:	100-300



