

## 兔抗 ASNS 多克隆抗体

中文名称：兔抗 ASNS 多克隆抗体

英文名称： Anti-ASNS rabbit polyclonal antibody

别名： TS11

相关类别： 一抗

抗原： ASNS

储存： 冷冻（-20℃）

宿主： Rabbit

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	The protein encoded by this gene is involved in the synthesis of asparagine. This gene complements a mutation in the temperature-sensitive hamster mutant ts11, which blocks progression through the G1 phase of the cell cycle at nonpermissive temperature. Alternatively spliced transcript variants have been described for this gene.
<b>Applications:</b>	ELISA, WB, IHC
<b>Name of antibody:</b>	ASNS
<b>Immunogen:</b>	Fusion protein of human ASNS
<b>Full name:</b>	asparagine synthetase (glutamine-hydrolyzing)
<b>Synonyms :</b>	TS11

<b>SwissProt:</b>	P08243
<b>ELISA Recommended dilution:</b>	1000-5000
<b>IHC positive control:</b>	Human colon cancer
<b>IHC Recommend dilution:</b>	50-200
<b>WB Predicted band size:</b>	64 kDa
<b>WB Positive control:</b>	Lovo and 231 cells
<b>WB Recommended dilution:</b>	500-2000



