

COX19 抗原（重组蛋白）

中文名称：COX19 抗原（重组蛋白）

英文名称：COX19 Antigen (Recombinant Protein)

别名：COX19 cytochrome c oxidase assembly factor; MGC104475

储存：冷冻（-20℃）

相关类别：抗原

概述

Full length fusion protein

技术规格

Full name:	COX19 cytochrome c oxidase assembly factor
Synonyms:	MGC104475
Swissprot:	Q49B96
Gene Accession:	BC110420
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	COX19 encodes a cytochrome c oxidase (COX)-assembly protein. The <i>S. cerevisiae</i> Cox19 protein may play a role in metal transport to the mitochondrial intermembrane space and assembly of complex IV of the mitochondrial respiratory chain.