

兔抗 ASPA 多克隆抗体

中文名称：兔抗 ASPA 多克隆抗体

英文名称： Anti-ASPA rabbit polyclonal antibody

别名： ASP; ACY2

相关类别： 一抗

抗原： ASPA

储存： 冷冻（-20℃）

宿主： Rabbit

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene encodes an enzyme that catalyzes the conversion of N-acetyl_L-aspartic acid (NAA) to a spartate and acetate. NAA is abundant in the brain where hydrolysis by aspartoacylase is thought to help maintain white matter. This protein is an NAA scavenger in other tissues. Mutations in this gene cause Canavan disease. Alternatively spliced transcript variants have been found for this gene.
Applications:	ELISA, WB, IHC
Name of antibody:	ASPA

Immunogen:	Fusion protein of human ASPA
Full name:	aspartoacylase
Synonyms :	ASP; ACY2
SwissProt:	P45381
ELISA Recommended dilution:	2000-10000
IHC positive control:	Human liver cancer and Human colon cancer
IHC Recommend dilution:	50-200
WB Predicted band size:	36 kDa
WB Positive control:	Mouse brain tissue
WB Recommended dilution:	1000-5000



