

CPNE3 抗原（重组蛋白）

中文名称： CPNE3 抗原（重组蛋白）

英文名称： CPNE3 Antigen (Recombinant Protein)

别名： copine 3; CPN3; PRO1071

储存： 冷冻（-20℃）

相关类别： 抗原

概述

Fusion protein corresponding to a region derived from 338-537 amino acids of human CPNE3

技术规格

Full name:	copine 3
Synonyms:	CPN3; PRO1071
Swissprot:	O75131
Gene Accession:	BC036242
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	Calcium-dependent membrane-binding proteins may regulate molecular events at the interface of the cell membrane and cytoplasm. This gene encodes a protein which contains two type II C2 domains in the amino-terminus and an A domain-like sequence in the carboxy-terminus. The A domain mediates interactions between integrins and extracellular ligands.