

兔抗 ASTN2 多克隆抗体

中文名称：兔抗 ASTN2 多克隆抗体

英文名称：Anti-ASTN2 rabbit polyclonal antibody

别名：astrotactin 2; bA67K19.1

抗原：ASTN2

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human, Mouse

相关类别：一抗

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes a protein that is expressed in the brain and may function in neuronal migration, based on functional studies of the related astrotactin 1 gene in human and mouse. A deletion at this locus has been associated with schizophrenia. Multiple transcript variants encoding different proteins have been found for this locus.
Applications:	ELISA, IHC
Name of antibody:	ASTN2
Immunogen:	Fusion protein of human ASTN2
Full name:	astrotactin 2

Synonyms:	bA67K19.1
SwissProt:	O75129
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human colorectal cancer
IHC Recommend dilution:	50-300



