

兔抗 COX18 多克隆抗体

中文名称：兔抗 COX18 多克隆抗体

英文名称： Anti-COX18 rabbit polyclonal antibody

别名： COX18HS

相关类别： 一抗

抗原： COX18

储存： 冷冻（-20℃）

宿主： Rabbit

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene encodes a cytochrome c oxidase assembly protein. The encoded protein is essential for integral membrane protein insertion into the mitochondrial inner membrane. It is also required for cytochrome c oxidase assembly and activity. Alternative splicing results in multiple transcript variants.
Applications:	ELISA, IHC
Name of antibody:	COX18
Immunogen:	Synthetic peptide of human COX18
Full name:	cytochrome c oxidase assembly factor COX18
Synonyms:	COX18HS

SwissProt:	Q8N8Q8
ELISA Recommended dilution:	1000-2000
IHC positive control:	Human thyroid cancer
IHC Recommend dilution:	10-50

