

兔抗 CPA2 多克隆抗体

中文名称：兔抗 CPA2 多克隆抗体

英文名称： Anti-CPA2 rabbit polyclonal antibody

相关类别： 一抗

储 存： 冷冻（-20℃）

抗 原： CPA2

宿 主： Rabbit

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	Three different forms of human pancreatic procarboxypeptidase A have been isolated. The encoded protein represents the A2 form, which is a monomeric protein with different biochemical properties from the A1 and A3 forms. The A2 form of pancreatic procarboxypeptidase acts on aromatic C-terminal residues and is a secreted protein.
Applications:	ELISA, WB, IHC
Name of antibody:	CPA2
Immunogen:	Fusion protein of human CPA2
Full name:	carboxypeptidase A2 (pancreatic)
SwissProt:	P48052
ELISA Recommended dilution:	2000-5000

IHC positive control:	Human liver cancer and human breast cancer
IHC Recommend dilution:	50-200
WB Predicted band size:	47 kDa
WB Positive control:	Mouse pancreas tissue
WB Recommended dilution:	500-2000



