

兔抗 ASXL2 多克隆抗体

中文名称：兔抗 ASXL2 多克隆抗体

英文名称： Anti-ASXL2 rabbit polyclonal antibody

别名： ASXH2

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： ASXL2

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	ASXL2 is a human homolog of the Drosophila asx gene. Drosophila asx is an enhancer of trithorax (see MIM 159555) and polycomb (see MIM 610231) (ETP) gene that encodes a chromatin protein with dual functions in transcriptional activation and silencing (Kato and Katoh, 2003 [PubMed 12888926]).
Applications:	ELISA, IHC
Name of antibody:	ASXL2
Immunogen:	Fusion protein of human ASXL2
Full name:	additional sex combs like 2, transcriptional regulator
Synonyms:	ASXH2

SwissProt:	Q76L83
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human cervical cancer
IHC Recommend dilution:	25-100

