

兔抗 CPA5 多克隆抗体

中文名称：兔抗 CPA5 多克隆抗体

英文名称： Anti-CPA5 rabbit polyclonal antibody

储 存： 冷冻（-20℃） 避光

抗 原： CPA5

宿 主： Rabbit

反应种属 Human

相关类别 一抗

标记物 Unconjugate

克隆类型 Unconjugate

技术规格

Background:	Carboxypeptidases have functions ranging from digestion of food to selective biosynthesis of neuroendocrine peptides. Members of the A/B subfamily of carboxypeptidases, such as CPA5, contain an approximately 90-amino acid pro region that assists in the folding of the active carboxypeptidase domain. Cleavage of the pro region activates the enzyme (Wei et al., 2002 [PubMed 11836249]).
Applications:	WB
Name of antibody:	CPA5
Immunogen:	Synthesized peptide derived from internal of human CPA5.
Full name:	carboxypeptidase A5

Synonyms :	0
SwissProt:	Q8WXQ8
WB Predicted band size:	49 kDa
WB Positive control:	293 cells lysate
WB Recommended dilution:	500-3000

