

CRB3 抗原（重组蛋白）

中文名称：CRB3 抗原（重组蛋白）

英文名称：CRB3 Antigen (Recombinant Protein)

别名：crumbs 3, cell polarity complex component

储存：冷冻（-20℃）

相关类别：抗原

概述

Fusion protein corresponding to a region derived from 27-123AA amino acids of human CRB3

技术规格

Full name:	crumbs 3, cell polarity complex component
Swissprot:	Q9BUF7
Gene Accession:	BC018409
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene encodes a member of the Crumbs family of proteins. This gene is widely expressed in epithelial tissues where the encoded protein isoforms play various roles such as the control of cytokinesis and ciliogenesis or the formation of tight junctions. Alternative splicing results in multiple transcript variants encoding different isoforms.