

CRCP 抗原（重组蛋白）

中文名称：CRCP 抗原（重组蛋白）

英文名称：CRCP Antigen (Recombinant Protein)

别名：CGRP receptor component; RCP; RCP9; RPC9; CGRPRCP; CGRP-RCP

储存：冷冻（-20℃）

相关类别：抗原

概述

Fusion protein corresponding to a region derived from 1-115 amino acids of human CRCP

技术规格

Full name:	CGRP receptor component
Synonyms:	RCP; RCP9; RPC9; CGRPRCP; CGRP-RCP
Swissprot:	O75575
Gene Accession:	BC040107
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene encodes a membrane protein that functions as part of a receptor complex for a small neuropeptide that increases intracellular cAMP levels. Alternate transcriptional splice variants, encoding different isoforms, have been characterized.