

## CRYZL1 抗原（重组蛋白）

中文名称：CRYZL1 抗原（重组蛋白）

英文名称：CRYZL1 Antigen (Recombinant Protein)

别名：4P11; QOH-1

储存：冷冻（-20℃）

相关类别：抗原

概述

Full length fusion protein

技术规格

<b>Full name:</b>	crystallin, zeta (quinone reductase)-like 1
<b>Synonyms:</b>	4P11; QOH-1
<b>Swissprot:</b>	O95825
<b>Gene Accession:</b>	BC013155
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system:</b>	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	This gene encodes a protein that has sequence similarity to zeta crystallin, also known as quinone oxidoreductase. This zeta crystallin-like protein also contains an NAD(P)H binding site. Alternatively spliced transcript variants have been observed but their full-length nature has not been completely determined.