

兔抗 CPN1 多克隆抗体

中文名称：兔抗 CPN1 多克隆抗体

英文名称： Anti-CPN1 rabbit polyclonal antibody

别名： CPN; SCPN

抗原： CPN1

储存： 冷冻（-20℃） 避光

宿主： Rabbit

反应种属： Human Mouse

相关类别： 一抗

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	Carboxypeptidase N is a plasma metallo-protease that cleaves basic amino acids from the C terminal of peptides and proteins. The enzyme is important in the regulation of peptides like kinins and anaphylatoxins, and has also been known as kininase-1 and anaphylatoxin inactivator. This enzyme is a tetramer comprised of two identical regulatory subunits and two identical catalytic subunits; this gene encodes the catalytic subunit. Mutations in this gene can be associated with angioedema or chronic urticaria resulting from carboxypeptidase N deficiency.
Applications:	WB

Name of antibody:	CPN1
Immunogen:	Synthesized peptide derived from C-terminal of human CPN1.
Full name:	carboxypeptidase N, polypeptide 1
Synonyms :	CPN; SCPN
SwissProt:	P15169
WB Predicted band size:	52 kDa
WB Positive control:	RAW264.7 cells lysate
WB Recommended dilution:	500-3000

