

CTNNBIP1 抗原（重组蛋白）

中文名称：CTNNBIP1 抗原（重组蛋白）

英文名称：CTNNBIP1 Antigen (Recombinant Protein)

别名：ICAT

储存：冷冻（-20℃）

相关类别：抗原

概述

Full length fusion protein

技术规格

Full name:	catenin, beta interacting protein 1
Synonyms:	ICAT
Swissprot:	Q9NSA3
Gene Accession:	BC014300
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	The protein encoded by this gene binds CTNNB1 and prevents interaction between CTNNB1 and TCF family members. The encoded protein is a negative regulator of the Wnt signaling pathway. Two transcript variants encoding the same protein have been found for this gene.