

## CTRC 抗原（重组蛋白）

中文名称： CTRC 抗原（重组蛋白）

英文名称： CTRC Antigen (Recombinant Protein)

别名： chymotrypsin C; CLCR; ELA4

储存： 冷冻（-20℃）

相关类别： 抗原

概述

Fusion protein corresponding to a region derived from 30-268 amino acids of human CTRC

技术规格

<b>Full name:</b>	chymotrypsin C
<b>Synonyms:</b>	CLCR; ELA4
<b>Swissprot:</b>	Q99895
<b>Gene Accession:</b>	BC015118
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system:</b>	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	This gene encodes a member of the peptidase S1 family. The encoded protein is a serum calcium-decreasing factor that has chymotrypsin-like protease activity. Alternatively spliced transcript variants have been observed, but their full-length nature has not been determined.