

兔抗 CPNE4 多克隆抗体

- 中文名称：兔抗 CPNE4 多克隆抗体
- 英文名称：Anti-CPNE4 rabbit polyclonal antibody
- 别名：copine 4 CPN4; COPN4
- 抗原：CPNE4
- 储存：冷冻 (-20℃)
- 宿主：Rabbit
- 相关类别：一抗
- 反应种属：Human, Mouse
- 标记物：Unconjugate
- 克隆类型：rabbit polyclonal

技术规格

Background:	This gene belongs to the highly conserved copine family. It encodes a calcium-dependent, phospholipid-binding protein, which may be involved in membrane trafficking, mitogenesis and development. Alternative splicing results in multiple transcript variants.
Applications:	ELISA, IHC
Name of antibody:	CPNE4
Immunogen:	Fusion protein of human CPNE4
Full name:	copine 4
Synonyms:	CPN4; COPN4
SwissProt:	Q96A23

ELISA Recommended dilution:	5000-10000
IHC positive control:	Human brain
IHC Recommend dilution:	20-100

