

兔抗 ATF7 多克隆抗体

中文名称：兔抗 ATF7 多克隆抗体

英文名称： Anti-ATF7 rabbit polyclonal antibody

别名： ATFA

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： ATF7

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	Plays important functions in early cell signaling. Binds the cAMP response element (CRE) (consensus: 5'-GTGACGT [AG][AG]-3'), a sequence present in many viral and cellular promoters. Activator of the NF-ELAM1/delta-A site of the E-selectin promoter. Has no intrinsic transcriptional activity, but activates transcription on formation of JUN or FOS heterodimers. Also can bind TRE promoter sequences when heterodimerized with members of the JUN family.
Applications:	ELISA, IHC
Name of antibody:	ATF7

Immunogen:	Full length fusion protein
Full name:	activating transcription factor 7
Synonyms:	ATFA
SwissProt:	P17544
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human lung cancer and human colorectal cancer
IHC Recommend dilution:	30-150



