

兔抗 CPNE8 多克隆抗体

中文名称：兔抗 CPNE8 多克隆抗体

英文名称：Anti-CPNE8 rabbit polyclonal antibody

别名：copine 8

储存：冷冻（-20℃）

抗原：CPNE8

宿主：Rabbit

相关类别：一抗

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	Calcium-dependent membrane-binding proteins may regulate molecular events at the interface of the cell membrane and cytoplasm. This gene is one of several genes that encode a calcium-dependent protein containing two N-terminal type II C2 domains and an integrin A domain-like sequence in the C-terminus.
Applications:	ELISA, IHC
Name of antibody:	CPNE8
Immunogen:	Fusion protein of human CPNE8
Full name:	copine 8
SwissProt:	Q86YQ8

ELISA Recommended dilution:	5000-10000
IHC positive control:	Human liver cancer
IHC Recommend dilution:	25-100

