

## CYB5R4 抗原（重组蛋白）

中文名称：CYB5R4 抗原（重组蛋白）

英文名称：CYB5R4 Antigen (Recombinant Protein)

别名：cytochrome b5 reductase 4; NCB5OR; cb5/cb5R; dJ676J13.1

储存：冷冻（-20℃）

相关类别：抗原

概述

Fusion protein corresponding to a region derived from 322-521 amino acids of human CYB5R4

技术规格

<b>Full name:</b>	cytochrome b5 reductase 4
<b>Synonyms:</b>	NCB5OR; cb5/cb5R; dJ676J13.1
<b>Swissprot:</b>	Q7L1T6
<b>Gene Accession:</b>	BC025380
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system:</b>	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	NADH-cytochrome b5 reductases are involved in desaturation and elongation of fatty acids, cholesterol biosynthesis, drug metabolism, and, in erythrocyte, methemoglobin reduction.