

## 兔抗 CPT1C 多克隆抗体

中文名称：兔抗 CPT1C 多克隆抗体

英文名称： Anti-CPT1C rabbit polyclonal antibody

别名： CATL1; CPT1P; CPTIC; SPG73; CPT1-B; CPTI-B

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： CPT1C

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a member of the carnitine/choline acetyltransferase family. The encoded protein regulates the beta-oxidation and transport of long-chain fatty acids into mitochondria, and may play a role in the regulation of feeding behavior and whole-body energy homeostasis. Alternatively spliced transcript variants encoding multiple protein isoforms have been observed for this gene.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	CPT1C
<b>Immunogen:</b>	Synthetic peptide of human CPT1C

<b>Full name:</b>	carnitine palmitoyltransferase 1C
<b>Synonyms :</b>	CATL1; CPT1P; CPTIC; SPG73; CPT1-B; CPTI-B
<b>SwissProt:</b>	Q8TCG5
<b>ELISA Recommended dilution:</b>	2000-5000
<b>IHC positive control:</b>	Human tonsil and Human esophagus cancer
<b>IHC Recommend dilution:</b>	25-100

