

## 兔抗 ATG16L1 多克隆抗体

中文名称: 兔抗 ATG16L1 多克隆抗体

英文名称: Anti-ATG16L1 rabbit polyclonal antibody

别名: IBD10; WDR30; APG16L; ATG16A; ATG16L

抗原: ATG16L1

储存: 冷冻 (-20℃)

宿主: Rabbit

反应种属: Human, Mouse

相关类别: 一抗

标记物: Unconjugate

克隆类型: rabbit polyclonal

### 技术规格

<b>Background:</b>	The protein encoded by this gene is part of a large protein complex that is necessary for autophagy, the major process by which intracellular components are targeted to lysosomes for degradation. Defects in this gene are a cause of susceptibility to inflammatory bowel disease type 10 (IBD10). Several transcript variants encoding different isoforms have been found for this gene.
<b>Applications:</b>	ELISA, WB, IHC
<b>Name of antibody:</b>	ATG16L1
<b>Immunogen:</b>	Synthetic peptide of human ATG16L1
<b>Full name:</b>	autophagy related 16 like 1

<b>Synonyms:</b>	IBD10; WDR30; APG16L; ATG16A; ATG16L
<b>SwissProt:</b>	Q676U5
<b>ELISA Recommended dilution:</b>	2000-5000
<b>IHC positive control:</b>	Human breast cancer and Human thyroid cancer
<b>IHC Recommend dilution:</b>	30-150
<b>WB Predicted band size:</b>	68 kDa
<b>WB Positive control:</b>	Hela and raji cell
<b>WB Recommended dilution:</b>	500-2000



