

兔抗 ATG16L2 多克隆抗体

中文名称：兔抗 ATG16L2 多克隆抗体

英文名称：Anti-ATG16L2 rabbit polyclonal antibody

别名：WDR80; ATG16B

相关类别：一抗

储存：冷冻（-20℃）避光

宿主：Rabbit

抗原：ATG16L2

反应种属：Human Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	May play a role in autophagy.
Applications:	WB
Name of antibody:	ATG16L2
Immunogen:	Synthesized peptide derived from internal of human ATG16L2.
Full name:	autophagy related 16 like 2
Synonyms:	WDR80; ATG16B
SwissProt:	Q8NAA4
WB Predicted band size:	69 kDa
WB Positive control:	NIH/3T3 cells and Mouse brain tissue lysates
WB Recommended dilution:	500-3000

