

## CYP4F12 抗原(重组蛋白)

中文名称: CYP4F12 抗原 (重组蛋白)

英文名称: CYP4F12 Antigen (Recombinant Protein)

别 名: cytochrome P450 family 4 subfamily F member 12; CYPIVF12; F22329\_1

储存: 冷冻 (-20℃)

相关类别: 抗原

概述

Fusion protein corresponding to a region derived from 325-524 amino acids of human CYP4F12

## 技术规格

Full name:	cytochrome P450 family 4 subfamily F member 12
Synonyms:	CYPIVF12; F22329_1
Swissprot:	Q9HCS2
Gene Accession:	BC035350
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene encodes a member of the cytochrome P450 superfamily of enzymes. The cytochrome P450 proteins are monooxygenases which ca talyze many reactions involved in drug metabolism and synthesis of ch olesterol, steroids and other lipids. This protein likely localizes to the e ndoplasmic reticulum. When expressed in yeast the enzyme is capable of oxdizing arachidonic acid. It can also catalyze the epoxidation of 22 :6n-3 and 22:5n-3 polyunsaturated long-chain fatty acids. This gene is



part of a cluster of cytochrome P450 genes on chromosome 19. Alter
native splicing results in multiple transcript variants.