

## 兔抗 CRHR1 多克隆抗体

中文名称：兔抗 CRHR1 多克隆抗体

英文名称： Anti-CRHR1 rabbit polyclonal antibody

别名： CRF1; CRHR; CRF-R; CRFR1; CRF-R1; CRFR-1; CRH-R1; CRHR1L; CRHR1f; CRF-R-1; CRH-R-1; CRH-R1h

相关类别： 一抗

抗原： CRHR1

储存： 冷冻（-20℃）

宿主： Rabbit

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

**Background:**

This gene encodes a G-protein coupled receptor that binds neuropeptides of the corticotropin releasing hormone family that are major regulators of the hypothalamic-pituitary-adrenal pathway. The encoded protein is essential for the activation of signal transduction pathways that regulate diverse physiological processes including stress, reproduction, immune response and obesity. Alternative splicing results in multiple transcript variants, one of which represents a read-through transcript with the neighboring gene MGC57346.

<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	CRHR1
<b>Immunogen:</b>	Synthetic peptide of human CRHR1
<b>Full name:</b>	corticotropin releasing hormone receptor 1
<b>Synonyms :</b>	CRF1; CRHR; CRF-R; CRFR1; CRF-R1; CRFR-1; CRH-R1; CRHR1L; CRHR1f; CRF-R-1; CRH-R-1; CRH-R1h
<b>SwissProt:</b>	P34998
<b>ELISA Recommended dilution:</b>	1000-2000
<b>IHC positive control:</b>	Human thyroid cancer
<b>IHC Recommend dilution:</b>	15-50

