

CYP4X1 抗原（重组蛋白）

中文名称：CYP4X1 抗原（重组蛋白）

英文名称：CYP4X1 Antigen (Recombinant Protein)

别名：cytochrome P450 family 4 subfamily X member 1; CYP4X1

储存：冷冻（-20℃）

相关类别：抗原

概述

Fusion protein corresponding to a region derived from 330-509 amino acids of human CYP4X1

技术规格

Full name:	cytochrome P450 family 4 subfamily X member 1
Synonyms:	CYP4X1
Swissprot:	Q8N118
Gene Accession:	BC028102
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene encodes a member of the cytochrome P450 superfamily of enzymes and is located within a cluster of genes belonging to this superfamily on chromosome 1. The cytochrome P450 proteins are monoxygenases which catalyze many reactions involved in drug metabolism and synthesis of cholesterol, steroids and other lipids. The expression pattern of a similar rat protein suggests that this protein may be involved in neurovascular function in the brain. Alternative splicing res

ults in multiple transcript variants.