

## 兔抗 ATOX1 多克隆抗体

中文名称：兔抗 ATOX1 多克隆抗体

英文名称： Anti-ATOX1 rabbit polyclonal antibody

别名： ATX1; HAH1

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： ATOX1

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a copper chaperone that plays a role in copper homeostasis by binding and transporting cytosolic copper to ATPase proteins in the trans-Golgi network for later incorporation to the ceruloplasmin. This protein also functions as an antioxidant against superoxide and hydrogen peroxide, and therefore, may play a significant role in cancer carcinogenesis. Because of its cytogenetic location, this gene represents a candidate gene for 5q-syndrome.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	ATOX1

<b>Immunogen:</b>	Fusion protein of human ATOX1
<b>Full name:</b>	antioxidant 1 copper chaperone
<b>Synonyms :</b>	ATX1; HAH1
<b>SwissProt:</b>	O00244
<b>ELISA Recommended dilution:</b>	2000-5000
<b>IHC positive control:</b>	Human breast cancer and human lung cancer
<b>IHC Recommend dilution:</b>	25-100

