

兔抗 Atp1a1(Ab-23) 多克隆抗体

中文名称：兔抗 Atp1a1(Ab-23) 多克隆抗体

英文名称： Anti-Atp1a1(Ab-23) rabbit polyclonal antibody

别名： Nkaa1b

相关类别： 一抗

储存： 冷冻（-20℃） 避光

宿主： Rabbit

抗原： Atp1a1(Ab-23)

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This is the catalytic component of the active enzyme, which catalyzes the hydrolysis of ATP coupled with the exchange of sodium and potassium ions across the plasma membrane. This action creates the electrochemical gradient of sodium and potassium ions, providing the energy for active transport of various nutrients.
Applications:	WB, IF
Name of antibody:	Atp1a1(Ab-23)
Immunogen:	Synthesized non-phosphopeptide derived from human ATP 1 α 1/Na+K+ ATPase1 around the phosphorylation site of serine 23 (K-K-S(p)-K-K).

Full name:	ATPase, Na ⁺ /K ⁺ transporting, alpha 1 polypeptide
Synonyms :	Nkaa1b
SwissProt:	P06685
WB Predicted band size:	113 kDa
WB Positive control:	Hela cells lysate
WB Recommended dilution:	500-3000
IF positive control:	Hela cells
IF Recommend dilution:	100-500



