

兔抗 CRMP1 多克隆抗体

中文名称：兔抗 CRMP1 多克隆抗体

英文名称： Anti-CRMP1 rabbit polyclonal antibody

别名： DRP1, DRP-1, CRMP-1, DPYSL1, ULIP-3

抗原： CRMP1

储存： 冷冻（-20℃）

宿主： Rabbit

相关类别： 一抗

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene encodes a member of a family of cytosolic phosphoproteins expressed exclusively in the nervous system. The encoded protein is thought to be a part of the semaphorin signal transduction pathway implicated in semaphorin-induced growth cone collapse during neural development. Alternative splicing results in multiple transcript variants.
Applications:	ELISA, IHC
Name of antibody:	CRMP1
Immunogen:	Fusion protein of human CRMP1
Full name:	collapsin response mediator protein 1

Synonyms :	DRP1, DRP-1, CRMP-1, DPYSL1, ULIP-3
SwissProt:	Q14194
ELISA Recommended dilution:	1000-5000
IHC positive control:	Human ovarian cancer and human lung cancer
IHC Recommend dilution:	25-100

