

## 兔抗 CRP 多克隆抗体

中文名称：兔抗 CRP 多克隆抗体

英文名称： Anti-CRP rabbit polyclonal antibody

别名： PTX1

抗原： CRP

储存： 冷冻（-20℃）

宿主： Rabbit

相关类别： 一抗

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	The protein encoded by this gene belongs to the pentaxin family. It is involved in several host defense related functions based on its ability to recognize foreign pathogens and damaged cells of the host and to initiate their elimination by interacting with humoral and cellular effector systems in the blood. Consequently, the level of this protein in plasma increases greatly during acute phase response to tissue injury, infection, or other inflammatory stimuli.
<b>Applications:</b>	ELISA, WB, IHC
<b>Name of antibody:</b>	CRP

<b>Immunogen:</b>	Fusion protein of human CRP
<b>Full name:</b>	C-reactive protein, pentraxin-related
<b>Synonyms :</b>	PTX1
<b>SwissProt:</b>	P02741
<b>ELISA Recommended dilution:</b>	2000-5000
<b>IHC positive control:</b>	Human thyroid cancer
<b>IHC Recommend dilution:</b>	50-200
<b>WB Predicted band size:</b>	25 kDa
<b>WB Positive control:</b>	Human liver cancer tissue
<b>WB Recommended dilution:</b>	500-2000



