

兔抗 CRTAM 多克隆抗体

中文名称：兔抗 CRTAM 多克隆抗体

英文名称： Anti-CRTAM rabbit polyclonal antibody

别名： CD355

抗原： CRTAM

储存： 冷冻（-20℃）

宿主： Rabbit

相关类别： 一抗

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	The CRTAM gene is upregulated in CD4 (see MIM 186940)-positive and CD8 (see CD8A; MIM 186910)-positive T cells and encodes a type I transmembrane protein with V and C1-like Ig domains. Interaction with CADM1 promotes natural killer (NK) cell cytotoxicity and interferon-gamma (IFN-gamma) secretion by CD8+ cells in vitro as well as NK cell-mediated rejection of tumors expressing CADM3 in vivo.
Applications:	ELISA, IHC
Name of antibody:	CRTAM
Immunogen:	Fusion protein of human CRTAM

Full name:	cytotoxic and regulatory T cell molecule
Synonyms :	CD355
SwissProt:	O95727
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human thyroid cancer
IHC Recommend dilution:	25-100

