

## 兔抗 CRTC2 多克隆抗体

中文名称：兔抗 CRTC2 多克隆抗体

英文名称：Anti-CRTC2 rabbit polyclonal antibody

别名：TORC2; TORC-2

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：CRTC2

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a member of the transducers of regulated cAMP response element-binding protein activity family of transcription coactivators. These proteins promote the transcription of genes targeted by the cAMP response element-binding protein, and therefore play an important role in many cellular processes. Under basal conditions the encoded protein is phosphorylated by AMP-activated protein kinase or the salt-inducible kinases and is sequestered in the cytoplasm.
<b>Applications:</b>	ELISA, WB

<b>Name of antibody:</b>	CRTC2
<b>Immunogen:</b>	Synthetic peptide of human CRTC2
<b>Full name:</b>	CREB regulated transcription coactivator 2
<b>Synonyms :</b>	TORC2; TORC-2
<b>SwissProt:</b>	Q53ET0
<b>ELISA Recommended dilution:</b>	1000-2000
<b>WB Predicted band size:</b>	73 kDa
<b>WB Positive control:</b>	Raji cells
<b>WB Recommended dilution:</b>	200-1000

