

兔抗 CRYAA 多克隆抗体

中文名称：兔抗 CRYAA 多克隆抗体

英文名称： Anti-CRYAA rabbit polyclonal antibody

别名： alpha-crystallin A chain

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： CRYAA

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	Contributes to the transparency and refractive index of the lens. Has chaperone-like activity, preventing aggregation of various proteins under a wide range of stress conditions.
Applications:	ELISA, WB
Name of antibody:	CRYAA
Immunogen:	Fusion protein of human CRYAA
Full name:	alpha-crystallin A chain
SwissProt:	P02489
ELISA Recommended dilution:	5000-10000
WB Predicted band size:	20 kDa

WB Positive control:	Human kidney tissue lysate
WB Recommended dilution:	200-1000

