

DCTN5 抗原（重组蛋白）

中文名称：DCTN5 抗原（重组蛋白）

英文名称：DCTN5 Antigen (Recombinant Protein)

储存 冷冻（-20℃）

相关类别 抗原

概述

Full length fusion protein

技术规格

Full name:	dynactin 5 (p25)
Swissprot:	Q9BTE1
Gene Accession:	BC004191
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene encodes a subunit of dynactin, a component of the cytoplasmic dynein motor machinery involved in minus-end-directed transport. The encoded protein is a component of the pointed-end subcomplex and is thought to bind membranous cargo. A pseudogene of this gene is located on the long arm of chromosome 1. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene.