

兔抗 ATP5PO 多克隆抗体

中文名称：兔抗 ATP5PO 多克隆抗体

英文名称： Anti-ATP5PO rabbit polyclonal antibody

别名： ATP synthase peripheral stalk subunit OSCP; ATPO;

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： ATP5PO

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	The protein encoded by this gene is a component of the F-type ATPase found in the mitochondrial matrix. F-type ATPases are composed of a catalytic core and a membrane proton channel. The encoded protein appears to be part of the connector linking these two components and may be involved in transmission of conformational changes or proton conductance. [provided by RefSeq, Jul 2008]
Applications:	ELISA, WB, IHC
Name of antibody:	ATP5PO
Immunogen:	Fusion protein of human ATP5PO
Full name:	ATP synthase peripheral stalk subunit OSCP

Synonyms:	ATPO; OSCP; ATP5O; HMC08D05
SwissProt:	P48047
IHC positive control:	Human thyroid cancer
IHC Recommend dilution:	100-300
WB Predicted band size:	23 kDa
WB Positive control:	SKOV3, Hela cell, Mouse liver tissue, Mouse heart tissue, Human placenta tissue, HepG2 cell, Human heart tissue lysates
WB Recommended dilution:	1000-5000



