

## 兔抗 CSF1 多克隆抗体

中文名称：兔抗 CSF1 多克隆抗体

英文名称： Anti-CSF1 rabbit polyclonal antibody

别名： MCSF; CSF-1

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： CSF1

反应种属： Human, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	The protein encoded by this gene is a cytokine that controls the production, differentiation, and function of macrophages. The active form of the protein is found extracellularly as a disulfide-linked homodimer, and is thought to be produced by proteolytic cleavage of membrane-bound precursors. The encoded protein may be involved in development of the placenta. Alternate splicing results in multiple transcript variants.
<b>Applications:</b>	ELISA, WB
<b>Name of antibody:</b>	CSF1
<b>Immunogen:</b>	Synthetic peptide of human CSF1

<b>Full name:</b>	colony stimulating factor 1 (macrophage)
<b>Synonyms:</b>	MCSF; CSF-1
<b>SwissProt:</b>	P09603
<b>ELISA Recommended dilution:</b>	5000-10000
<b>WB Predicted band size:</b>	60 kDa
<b>WB Positive control:</b>	Human heart tissue lysate
<b>WB Recommended dilution:</b>	200-1000

