

## 兔抗 CSF2RA 多克隆抗体

中文名称：兔抗 CSF2RA 多克隆抗体

英文名称： Anti-CSF2RA rabbit polyclonal antibody

别名： GMR; CD116; CSF2R; SMDP4; CDw116; CSF2RX; CSF2RY; GMCSFR; CSF2RAX; CSF2RAY; GM-CSF-R-alpha

抗原： CSF2RA

储存： 冷冻 (-20℃)

宿主： Rabbit

相关类别： 一抗

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

**Background:**

The protein encoded by this gene is the alpha subunit of the heterodimeric receptor for colony stimulating factor 2, a cytokine which controls the production, differentiation, and function of granulocytes and macrophages. The encoded protein is a member of the cytokine family of receptors. This gene is found in the pseudoautosomal region (PAR) of the X and Y chromosomes. Multiple transcript variants encoding different isoforms have been found for this gene, with some of the isoforms being membrane-bound and others being soluble.

<b>Applications:</b>	ELISA, WB, IHC
<b>Name of antibody:</b>	CSF2RA
<b>Immunogen:</b>	Synthetic peptide of human CSF2RA
<b>Full name:</b>	colony stimulating factor 2 receptor, alpha, low-affinity
<b>Synonyms :</b>	GMR; CD116; CSF2R; SMDP4; CDw116; CSF2RX; CSF2RY; GMCSFR; CSF2RAX; CSF2RAY; GM-CSF-R-alpha
<b>SwissProt:</b>	P15509
<b>ELISA Recommended dilution:</b>	2000-5000
<b>IHC positive control:</b>	Human cervical cancer and Human thyroid cancer
<b>IHC Recommend dilution:</b>	50-200
<b>WB Predicted band size:</b>	46 kDa
<b>WB Positive control:</b>	Human placenta tissue and HepG2 cells
<b>WB Recommended dilution:</b>	500-2000



