

## 兔抗 CSF2RB 多克隆抗体

中文名称：兔抗 CSF2RB 多克隆抗体

英文名称：Anti-CSF2RB rabbit polyclonal antibody

别名：CD131; IL3RB; IL5RB; SMDP5; CDw131

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：CSF2RB

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	The protein encoded by this gene is the common beta chain of the high affinity receptor for IL-3, IL-5 and CS F. Defects in this gene have been reported to be associated with protein alveolar proteinosis (PAP).
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	CSF2RB
<b>Immunogen:</b>	Synthetic peptide of human CSF2RB
<b>Full name:</b>	colony stimulating factor 2 receptor, beta, low-affinity (granulocyte-macrophage)
<b>Synonyms :</b>	CD131; IL3RB; IL5RB; SMDP5; CDw131
<b>SwissProt:</b>	P32927

<b>ELISA Recommended dilution:</b>	1000-2000
<b>IHC positive control:</b>	Human thyroid cancer and Human liver cancer
<b>IHC Recommend dilution:</b>	25-100

