

DHRS7 抗原（重组蛋白）

中文名称：DHRS7 抗原（重组蛋白）

英文名称：DHRS7 Antigen (Recombinant Protein)

别名：dehydrogenase/reductase 7; CGI-86; SDR34C1; retDSR4; retSDR4

储存：冷冻（-20℃）

相关类别：抗原

概述

Fusion protein corresponding to a region derived from 60-339 amino acids of human DHRS7

技术规格

Full name:	dehydrogenase/reductase 7
Synonyms:	CGI-86; SDR34C1; retDSR4; retSDR4
Swissprot:	Q9Y394
Gene Accession:	BC000637
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene encodes a member of the short-chain dehydrogenases/reductases (SDR) family, which has over 46,000 members. Members in this family are enzymes that metabolize many different compounds, such as steroid hormones, prostaglandins, retinoids, lipids and xenobiotics.