

## DIMT1 抗原（重组蛋白）

中文名称：DIMT1 抗原（重组蛋白）

英文名称： DIMT1 Antigen (Recombinant Protein)

别名： DIM1 dimethyladenosine transferase 1 homolog; DIM1; DIMT1L; HUSSY5; HSA9761

储存： 冷冻（-20℃）

相关类别： 抗原

概述

Fusion protein corresponding to a region derived from 1-313 amino acids of human DIMT1

技术规格

<b>Full name:</b>	DIM1 dimethyladenosine transferase 1 homolog
<b>Synonyms:</b>	DIM1; DIMT1L; HUSSY5; HSA9761
<b>Swissprot:</b>	Q9UNQ2
<b>Gene Accession:</b>	BC010874
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system:</b>	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	Specifically dimethylates two adjacent adenosines in the loop of a conserved hairpin near the 3'-end of 18S rRNA in the 40S particle.