

兔抗 CSF3 多克隆抗体

中文名称: 兔抗 CSF3 多克隆抗体

英文名称: Anti-CSF3 rabbit polyclonal antibody

别名: GCSF, CSF3OS, C17orf33

相关类别: 一抗

储存: 冷冻 (-20℃)

宿主: Rabbit

抗原: CSF3

反应种属: Human

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	The protein encoded by this gene is a cytokine that controls the production, differentiation, and function of granulocytes. The active protein is found extracellularly. Alternatively spliced transcript variants have been described for this gene.
Applications:	ELISA, IHC
Name of antibody:	CSF3
Immunogen:	Synthetic peptide of human CSF3
Full name:	colony stimulating factor 3 (granulocyte)
Synonyms :	GCSF, CSF3OS, C17orf33
SwissProt:	P09919

ELISA Recommended dilution:	1000-2000
IHC positive control:	Human esophagus cancer and Human tonsil
IHC Recommend dilution:	15-50

