

## DNAJB4 抗原（重组蛋白）

中文名称：DNAJB4 抗原（重组蛋白）

英文名称：DNAJB4 Antigen (Recombinant Protein)

别名：DjB4; HLJ1; DNAJW

储存：冷冻（-20℃）

相关类别：抗原

概述

Full length fusion protein

技术规格

<b>Full name:</b>	DnaJ (Hsp40) homolog, subfamily B, member 4
<b>Synonyms:</b>	DjB4; HLJ1; DNAJW
<b>Swissprot:</b>	Q9UDY4
<b>Gene Accession:</b>	BC034721
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system:</b>	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	DNAJB4 belongs to the evolutionarily conserved DNAJ/HSP40 protein family. For background information on the DNAJ family, see MIM 608375.